

Suro pedial

Company

Identification

Sex Age Weight Height

Right side Left side

Measurement date

Desired deadline

Applicator





Elements provided

- Measurements form
- Digital file
- Moulding with measurements form
- Profile/front photos with measurements form


Information

- Foot riser
- Articulated boot
- Rigid boot
- Night
- Day

Foot angulation

- Varus 
- Valgus 
- Dorsiflexion 
- Plantar flexion 

Remarks



Box reserved for CSM

Apparatus No
Operator

CSM production ref.

Suro pedial

 Date

 Request for quote

 Order

 Company

 Delivery address

 Applicator

 Private label brand

 Desired deadline

 Remarks
Elements to be provided for modelling

Moulding with measurements form
OR
 Measurements form
OR
 Profile/front photos with measurements form
OR
 SL digital file (Stereo-Lithography)

Elements to be provided for milling

Rodin4D digital file
OR
 SL digital file (Stereo-Lithography)
OR
 Captevia photo (.cpv)

Material	Thickness
HP polyethylene	
Polypropylene	
Polypropylene copolymer	

Padded foam	Colour	Thickness
Plastazote		
3D fabric		

Transfer paper	
Indicate the desired reference <i>(Refer to our Supplies catalogue or ask for our colour-sampler document)</i>	
<input type="checkbox"/> Transfer paper supplied (tick the box)	

Joint	
<input type="checkbox"/> Ankle	
<input type="checkbox"/> Other	

Velcro strap	Material	Colour
<input type="checkbox"/> PVC		
<input type="checkbox"/> PU		

References	Description	Order
SBS100	Mould scan	
SBR100	Boot grinding	
SBU100	Boot milling max H. 64 cm	
SOIMP	Preparation of a plastered counter-mould	

Composite		
SOI59CRL3	JAD blade carbon riser in 3 parts	
SOI59CRL	JAD blade carbon riser	
SOI59CAEL	JAD blade carbon anti-collapse	
SOI59CAR	Carbon boot, comprising 1 return joint	
SOI59LCC	Made-to-measure carbon composite JAD blade	
SOI59C	Carbon boot	

Polyolefin		
SOI59T	Thermoformed riser boot	
SOIACS	Flexible ankle joint on thermoformed boot	
SOI59ODAP1	Foot Absolute Discharge orthosis	
SOI59ODAP2	Foot Partial Discharge orthosis	

ORTHOSIS FINISHING		
SOI59G	Fabric trimming + velcro	

Box reserved for CSM

 Apparatus No
 Operator
 Delivery date

 Remarks